

## “D” IMPLANT DIMENSIONS

Implant	Mesial-distal width at crest	Neck Diameter	Buccal-lingual width at crest	Length
D1 short	4.0 mm	4.0 mm	3.5 mm	16 mm
D2 short	5.5 mm	4.5 mm	3.5 mm	16 mm
D3 short	7.5 mm	4.5 mm	3.5 mm	16 mm
D4 short	10 mm	4.5 mm	3.5 mm	16 mm
D1	4.0 mm	4.0 mm	3.5 mm	20 mm
D2	5.5 mm	4.5 mm	3.5 mm	20 mm
D3	7.5 mm	4.5 mm	3.5 mm	20 mm
D4	10 mm	4.5 mm	3.5 mm	20 mm