

Traus SUS10



Ultrasonic + Implant Dual Engine Stable and Powerful Output

- Offers an Ultrasonic Piezo & Implant motor in one unit
- Compact size
- Affordable and competitively price
- Easy to read control panel with vivid LCD
- Optic or non-optic motor are available
- Includes contra angle handpieces & 10 Piezo tips

- Complete set includes controller and pedal, 1 piezo handpiece, 1 motor, one 20:1 contra angle handpiece, inner irrigation clip, 2 handpiece stands, 10 Piezo tips with tip holder, handpiece cage, irrigation tube, tube holders, 2 poles for irrigation and accessories
- Motor auto-calibration function
- E-type motor with ISO standard connection
- Optional 32:1 angle handpiece
- Maximum speed at 2,000rpm with 20:1, 1250rpm with 32:1 angle handpiece
- Max 70Ncm with 20:1 handpiece, 80Ncm with 32:1 handpiece
- Actual speed and torque display
- Automatic overload protection system
- 9 programs memory function
- Motor & angle are autoclavable



Effectiveness of Piezo handpiece

Tip Scanning Function

bars on display to show the tip usage and tip connection condition to handpiece

Operation Standby Mode

piezo handpiece operating on loaded condition which extends the span life of Piezo handpiece and tip

Optimized Handpiece Torque Efficiency

checks the condition of the handpiece & connected tip to optimize output

Traus SIP10



Optic

Optic or non-optic motor & handpiece available

- Complete set includes controller and pedal, one 40,000rpm BLDC motor, one 20:1 angle handpiece, one handpiece stand, inner irrigation clip, tube holders, irrigation pole
- Motor auto-calibration function
- E-type motor with ISO standard connection
- Optional 32:1 angle handpiece available
- Maximum speed at 2,000rpm with 20:1, 1250rpm with 32:1 angle handpiece
- Max 70Ncm with 20:1 handpiece, 80Ncm with 32:1 handpiece
- Actual speed and torque displayed
- Automatic overload protection system
- 9 programs memory function
- Motor & angle are autoclavable

Implant CUBE



Economical and simple to use implant motor

- Complete set includes controller and pedal, one 50,000rpm BLDC motor, one 20:1 angle handpiece, one handpiece stand, inner irrigation clip, tube holders, irrigation pole
- Optional 32:1 angle handpiece available
- Maximum speed at 2,500rpm with 20:1, 1560rpm with 32:1 angle handpiece
- Max 55Ncm with 20:1 handpiece, 65Ncm with 32:1 handpiece
- Actual speed & torque displayed
- 9 programs memory function
- Motor & angle are autoclavable



TRAUS ENDO



www.saeshin.com

SAESHIN AMERICA, INC.
216 Technology Dr . Ste F, Irvine, CA 92618

T. 949-825-6925 F. 949-825-6926
E. saeshinamerica@saeshin.com